
Vacotron S

Bedieningshandleiding
Operating Instructions
Gebrauchsanweisung
Mode d'emploi
Instrucciones de manejo

Copyright:

Enraf-Nonius B.V.
P.O. Box 12080
3004 GB ROTTERDAM
The Netherlands
Tel: +31 (0)10 – 20 30 600
Fax: +31 (0)10 – 20 30 699
info@enraf-nonius.nl
www.enraf-nonius.com

Part number: 1452.750-41
2nd edition, February 2006

Vacotron S

Bedieningshandleiding
Operating Instructions
Gebrauchsanweisung
Mode d'emploi
Instrucciones de manejo

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Table of Contents

1	Foreword	14
2	Introduction	14
	2.1 Installation	14
	2.2 Technical maintenance	15
	2.3 Product Liability	15
3	Starting to use the unit	16
	3.1 Connecting the vacuum hoses (left side)	16
	3.2 Patient cables (left side)	16
	3.3 Mains adapter (backside).....	16
4	Controls	16
5	Operating sequence	18
	5.1 Vacuum electrodes.....	18
	5.1.1 <i>Switching on</i>	18
	5.1.2 <i>Positioning the electrodes</i>	18
	5.1.3 <i>Correcting the vacuum</i>	18
	5.1.4 <i>Selecting the current</i>	18
	5.1.5 <i>Ending treatment</i>	18
	5.2 Flexible rubber electrodes	18
6	Contra-Indications	19
7	Maintenance by the User	19
8	Fault conditions	20
9	Specifications and Technical data	20
10	Ordering data	20
	10.1 Standard accessories	20
	10.2 Additional accessories	20

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1 Foreword

The Vacotron S pulsed vacuum unit provides fast and simple electrode positioning using vacuum electrodes. Furthermore, it is extremely simple to switch between treatment with vacuum electrodes and treatment using flexible rubber electrodes. As a result, the Vacotron S is an ideal aid to electrotherapy.

When vacuum electrodes are applied, effective use is made of the whole electrode area, as the electrodes make good contact with the skin. The massage effect resulting from the pulsed vacuum ensures a good flow of blood through the skin under the electrode. This reduces the resistance of the skin, increases the effectiveness of the stimulating current, and reduces the treatment time.

Constant or pulsed vacuum

When a constant vacuum strength is used, the vacuum pulse switches off as soon as the selected vacuum strength is reached, making the Vacotron S virtually silent in operation. A control system ensures that the vacuum remains constant. The pump switches on for short intervals during the treatment to maintain the vacuum.

When pulsed vacuum is used, the pump is switched on and off intermittently. The pulsating rhythm of the vacuum massage is continuously or pulsating (two possibilities) adjustable. Pulsating vacuum is very useful for patients who are sensitive to the electric current, as the sensation of the current tends to be suppressed.

Electrodes

There is a choice of large and small electrodes. The areas of the electrodes correspond to those of the 4 x 6 cm and 6 x 8 cm flexible rubber electrodes. The vacuum electrodes are sufficiently flexible to ensure optimum contact with the skin, but rigid enough to prevent any changes in the contour of the part being treated, allowing full advantage to be taken of the massage effect of the pulsed vacuum.

Smart connections for simple operation

The Vacotron S has been designed to ensure easy operation. This is reflected in the design of the electrode connections. The Vacotron S has connections for both vacuum electrodes and flexible rubber electrodes, providing a permanently available choice between the two types without the need to change the cables. The vacuum electrodes are automatically selected as soon as the vacuum is set.

Finally

The Vacotron S has been designed and manufactured to ensure the highest level of safety, and the unit fully complies with the stringent IEC 60601-1 standards. We are confident that your unit will continue to give satisfaction over many years of use. Nevertheless, if you have any queries or suggestions, please contact your authorized distributor.

2 Introduction

It is important that you read these operating instructions carefully before using the Vacotron S. Please make sure that these instructions are available to all personnel who operate the equipment.

Pay attention to the following before using the Vacotron S:

1. Keep yourself informed of the contra-indications (see further).
2. The apparatus may not be used in close proximity (i.e. less than 2 metres) to shortwave or microwave equipment.

The manufacturer cannot be held responsible for the results of using this apparatus for any purposes other than described in these operating instructions.

2.1 Installation

- Do not install the unit in a location near a heat source such as a radiator.
- Avoid exposure to direct sunlight, rain, excessive dust, dampness, mechanical vibrations and shocks.
- This unit may not be used in so-called "wet rooms" (hydrotherapy rooms).

Should any liquid or foreign object get into the housing of this equipment, turn off the power switch immediately, unplug the unit from the wall socket and have it checked by an authorized expert.

Connection

- Connection to the electrical mains supply must comply with the national requirements regarding medical rooms.
- Before connection of this apparatus to the supply, check that the mains voltage and frequency stated on the type plate of this apparatus correspond to that of the mains supply.
- This apparatus may only be connected to an earthed wall socket.
- The unit should only be used with the mains cable supplied with this equipment.
- The mains supply is an integral part of the power circuit which is a fundamental part of the safety of this unit. The approval of the Vacotron S is only applicable in combination with adapter ENA-1550. It is not permitted to use any other adapter.

The Vacotron S can be used with a combination cable in connection to a Sonopuls 492 or Endomed 482e. In that case the adapter of these units are used.

2.2 Technical maintenance

We recommend to have the unit checked annually. This may be done by your supplier, or by an other agency, authorized by the manufacturer. It is also recommended that a record of the service history is kept for all activities relating to service and maintenance. In some countries this is even obligatory.

No attempt should be made to open the unit. Maintenance and all repairs should be carried out by an authorized agency.

The manufacturer will not be held responsible for the results of maintenance or repairs by unauthorized persons.

Electrolytic effects

Electrolysis occurs under the electrodes when Galvanic or currents with a DC component are applied. Because the greatest concentration of electrolytic by-products caused by ion migration occur under the electrodes, we recommend the use of the supplied sponges to keep the effects to a minimum. Make sure that the sponges are kept well moistened.

Current density

The particular standard for Therapeutic Electrical Stimulation Equipment is IEC 60601-2-10 (BS5724 part 2.10). The maximum permissible Current Density is 2 mA R.M.S. per cm².

To verify that this value is not exceeded in electrotherapy, divide the current output in mA by the effective electrode area in cm²; a value of less than two should result.

Note that placing the electrodes should be done carefully and solidly, and that the use of small electrodes in combination with high intensities may cause skin irritations or burns.

Electromagnetic interference

Simultaneous connection of a patient to High Frequency surgical equipment and an electrical stimulator must not be done as this may result in burns at the site of the stimulator electrodes.

Operation in close proximity (i.e. less than 2 metres) to working shortwave or microwave therapy unit may produce instability in the Vacotron S output.

To prevent Electromagnetic Interference, we strongly advise using separate mains groups (phases) for the Vacotron S and the shortwave or microwave equipment, and a distance of at least 2 metres is kept between the Vacotron S and the shortwave or microwave equipment. Ensure that the mains cable of the shortwave/microwave unit does not come near the Vacotron S or the patient.

If problems with Electromagnetic Interference persist, please contact your authorized distributor.

2.3 Product Liability

A law on Product Liability has become effective in many countries. This Product Liability law implies, amongst other things, that once a period of 10 years has elapsed after a product has been brought into circulation, the manufacturer can no longer be held responsible for possible shortages of the product.

To the maximum extent permitted by applicable law, in no event will Enraf-Nonius or its suppliers or resellers be liable for any indirect, special, incidental or consequential damages arising from the use of or inability to use the product, including, without limitation, damages for loss of goodwill, work and

productivity, computer failure or malfunction, or any and all other commercial damages or losses, even if advised of the possibility thereof, and regardless of the legal or equitable theory (contract, tort or otherwise) upon which the claim is based. In any case, Enraf-Nonius's entire liability under any provision of this agreement shall not exceed in the aggregate the sum of the fees paid for this product and fees for support of the product received by Enraf-Nonius under a separate support agreement (if any), with the exception of death or personal injury caused by the negligence of Enraf-Nonius to the extent applicable law prohibits the limitation of damages in such cases. Enraf-Nonius can not be held liable for any consequence resulting from incorrect information provided by it's personnel, or errors incorporated in this manual and / or other accompanying documentation (including commercial documentation) The opposing party (product's user or it's representative) shall disclaim Enraf-Nonius from all claims arising from third parties, whatever nature or whatever relationship to the opposing party.

3 Starting to use the unit

(For pictures and drawings please refer to page 47)

Connect the electrotherapy unit using cable (8) wich is located at the front panel of the Vacotron S.

3.1 Connecting the vacuum hoses (left side)

channel 1

- Connect the red hose to the left hand socket 2 (+).
- Connect the black hose to the right hand socket 2 (-).

channel 2

- Connect the red hose to the left hand socket 3 (+).
- Connect the black hose to the right hand socket 3 (-).

3.2 Patient cables (left side)

Connect the patient cables of your electrotherapy unit to connection 5 for channel 1 and connection 6 for channel 2.

3.3 Mains adapter (backside)

Connect the mains adapter of your therapy unit to channel 10. Connect the combinationcable 9 to your therapy unit

Attention: you can only use the connection cable in combination with a mains adapter type ENA-1550 (Endomed 482e and Sonopuls 492). In case other therpeutical devices are used this connection cable can not be used and the (optional) mains adapter ENA-1550 should be used in order to connect the Vacotron S

Attention: in case you use the Vacotron S in combintation with the Sonopuls 492, and you use both current channels of the Sonopuls 492 and the ultrasound channel at full power, than it is recommended not to use the connection cable. In that case we recommend to use the (optional) ENA-1550 adapter.

4 Controls

1 Vacuum regulator

This regulator sets the vacuum level, and is also used to switch the vacuum system on and off.

No vacuum:

Turn the regulator fully to the left.

The vacuum electrodes are inoperative. Indicator lamp 4 is out. The patient cables connected to connections 5 and 6 are operative.

Vacuum:

Turn the regulator to the right. The indicator lamp 4 by the vacuum electrodes is on, indicating that the vacuum electrodes are now the operative electrodes. Further rotation to the right increases the vacuum

2 Vacuum channel 1

Left hand connection: the red vacuum hose (+).

Right hand connection: the black vacuum hose (-).

Vacuum channel 1 is connected to one channel of the electrotherapy unit via connection 8.

3 Vacuum channel 2

Left hand connection: the red vacuum hose (+).

Right hand connection: the black vacuum hose (-).

Vacuum channel 2 is connected to the second channel of the electrotherapy unit via connection 8.

Connection of any cables other than the vacuum cables can adversely affect the safety of the patient and the functioning of the unit, and is therefore not permitted.

4 Vacuum electrodes indicator lamp

Lamp is lit when the vacuum electrodes are operative.

Lamp is off when the flexible rubber electrodes are operative.

5 Patient cable connection channel 1

Connect one of the patient cables supplied for use with the flexible rubber electrodes to this connection.

This connection is selected when no vacuum has been set. Channel 1 is connected to one channel of the electro-therapy unit via connection 8.

6 Patient cable connection channel 2

Connect one of the patient cables supplied for use with the flexible rubber electrodes to this connection.

This connection is selected when no vacuum has been set. Channel 2 is connected to one channel of the electro-therapy unit via connection 8.

Note: the use of other cables as the ones supplied with the unit can cause danger for the patient or influence proper functioning of the unit negatively. Therefore they are not permitted.

7 Key for constant vacuum, pulsed vacuum 1:1, pulsed vacuum 0,5:0,5

After switching the Vacotron to on it will supply constant vacuum (the constant vacuum indicator lamp lights up). After pressing the key 7 the Vacotron will supply pulsed vacuum with a fixed pulse:pause ratio of 1:1 second (the indicator lamp 1:1 lights up). After pressing key 7 again the Vacotron will supply pulsed vacuum with a fixed pulse:pause ratio of 0,5:0,5 seconds (the indicator lamp 0,5:0,5 lights up). Pressing key 7 again will result in constant vacuum (and the continuous vacuum lamp lights up).

8 Connection for the electrotherapy unit

This connection is used for connecting the patient-channel of the electrotherapy unit to the Vacotron S.

9 Connection cable for the power supply

This cable can be connected with the electrotherapy unit in case the adapter of this unit is a type ENA-1550.

10 Connection adaptor

This is the connection pin for the adaptor type ENA-1550

11 On-off switch

This switch is used to turn the Vacotron S on or off.

12 Connection tube for the water reservoir

The water reservoir ensures that any excess water from the sponges is collected, and cannot be blown out through the pump. Please refer to the section on Maintenance for emptying and cleaning the reservoir.

13 Closed connection tube for the water reservoir

This tube is connected to the other end of the connection tube for the water reservoir under normal circumstances.

14 Indication plate

Provides information on the unit, such as type and serial number, as well as connection data such as mains voltage and current consumption.

15 Warning sticker

This sticker indicates which adaptor must be used.

16 Air outlet

the air outlet is provided with a noise reduction filter.

5 Operating sequence

5.1 Vacuum electrodes

5.1.1 Switching on

- Switch on the Vacotron S and the electrotherapy unit.

5.1.2 Positioning the electrodes

- Choose the appropriate electrode size for the body part to be treated.
- Dampen the sponges. Squeeze them slightly to remove excess water, and place them in the vacuum electrodes.
- Select constant vacuum, this will ease the positioning of the electrodes.
- Place the vacuumelectrodes on a smooth surface (e.g. on a table).
- Turn the vacuumintensity up a little, so as to create a vacuum in the electrodes.
- Now move the electrodes one by one from the smooth surface to the body of the patient.
- Set the vacuum intensity to the desired value.
- Set the treatment time.
- The treatment starts.

5.1.3 Correcting the vacuum

- Correct the vacuum if necessary with regulator 1. A low vacuum is sufficient.
- If required, select the pulsed vacuum with the key 7.

5.1.4 Selecting the current

- Select the required current type, and gradually increase the current. Please refer to the Operating Instructions supplied with the electrotherapy unit.

5.1.5 Ending treatment

- Reduce the current of the electrotherapy unit to zero.
- Turn the vacuum regulator back to position 2, and remove the vacuum electrodes from the skin.
- Turn the vacuum regulator back to zero.

5.2 Flexible rubber electrodes

- The vacuum regulator should be in the off position. Indication light 4 is off
- Treatment continues as described in the Operating Instructions of the Electrotherapy unit.

N.B.

During treatment with the flexible rubber electrodes, simultaneous use of the vacuum electrodes is not possible. The vacuum regulator should remain in the zero position during treatment with the flexible rubber electrodes.

6 *Contra-Indications*

Supplementary pulsed vacuum therapy can be applied in almost every situation where low-frequency electrotherapy is indicated.

The use of the vacuum electrodes is only contra-indicated in the case of inflammatory processes, or where the risk of bleeding at the electrode sites cannot be excluded.

7 *Maintenance by the User*

Treatment unit

First switch off the unit and disconnect the mains adapter. Remove the connection cable of the electrotherapy unit. The unit can be cleaned with a damp cloth. Use lukewarm water and, if required, a liquid household detergent. Do not use abrasive cleaners.

Electrodes and sponges

The vacuum electrodes and sponges should be cleaned with lukewarm water. In the case of persistent dirt, and for disinfection, a 70% alcohol solution may be used.

Sponges should be replaced regularly. It is recommended to keep sponges and a spare electrode in stock.

Calcium scale can be deposited on the metal surfaces of the electrodes. This has an insulating effect. In order to maintain optimum conductivity, these surface should be regularly cleaned and polished.

Electrode cables

Check the cables regularly for damage and/or poor contacts. It is recommended to keep a spare electrode cable in stock.

Water reservoir

In case the water reservoir is full all lights will blink simultaneously and you will hear 5 beeps. The Vacotron S should then be switched off and the water reservoir should be emptied by means of disconnecting tube 12 from tube connector 13 and letting the liquid run out. After this the tube should be replaced and the unit can be turned on again.

Cleaning of the water reservoir and tubes.

Before you will run the cleaning procedure, the peripheral ends of the vacuum tubes should be placed in a cup with cleaning liquid (see hereunder) so the Vacotron S can be cleaned internally with this liquid. After this press button 7 for 3 seconds. The indication lights will now lit one by one and you will hear 3 beeps. After this the pump will automatically work in full power until the reservoir is full with cleaning liquid. In case the reservoir is full, than empty is as described above.

The following registered products may be used for disinfecting the water reservoir:

BAKTOLAN to 5%
CHINOSOL to 1%
CHLORAMIN solution
ELMOCID Gamma to 2%
MEFAROL to 1%.
MERCKOJOD to 1%
MERFEN
PERHYDROL
PERODIN
SAGROTAN to 2%
ZEPHIROL to 5%.

8 Fault conditions

Lamps fail to light

Check the mains adapter and the mains cables.

No patient current

Check that the correct electrodes have been selected. If the vacuum regulator is in the zero position, the flexible rubber electrodes are selected automatically. As soon as a vacuum is set, the vacuum electrodes are selected.

Check operation of the electrotherapy unit using the operating instructions provided (if necessary, with the Vacotron S disconnected).

Check the cables from the electrotherapy unit to the Vacotron S.

Patient current too low

The electrical resistance of the electrodes is too high. This can be due to the sponges not being damp enough, or to a layer of calcium scale or corrosion on the electrodes.

9 Specifications and Technical data

Vacuum	:	0-600 mbar, continuously adjustable
Pulsed vacuum	:	pulse : pause = 0,5:0,5 seconds, 1:1 second
No. of vacuum electrodes	:	max. 4
Safety class	:	I* type BF**, according to IEC 60601-1
Dimensions	:	24,5 x 21 x 9 cm
Weight	:	2,5 kg
Safety	:	TÜV Rheinland
		This device is in conformity with the European Directive (93/42/EEC).
Adaptor		
Type Adaptor	:	ENA-1550
Mains voltage	:	100 - 240 Volt $\pm$ 10%
Frequency	:	50/60 Hz.
DC Output	:	15 Volt
Max. Output current	:	3,3 A

*I: The equipment has a safety earth (ground) connection, and must be connected to an earthed (grounded) wall socket.

**BF: The equipment has a floating patient circuit. The Vacotron S meets the requirements of IEC 60601-1 and IEC 60601-2-10.

Technical modifications reserved.

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