
Endomed 182

Bedieningshandleiding
Operating Instructions
Gebrauchsanweisung
Mode d'emploi
Instrucciones de manejo



Copyright:

**Enraf
Nonius**

Enraf-Nonius B.V.
P.O. Box 12080
3004 GB Rotterdam
The Netherlands
Tel.: +31(0)10 - 20 30 600
Fax: +31(0)10 - 20 30 699

Article number : 1416.751-41
August 2005

Endomed 182

Bedieningshandleiding
Operating Instructions
Gebrauchsanweisung
Mode d'emploi
Instrucciones de manejo

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TABEL OF CONTENTS

Chapter 1 - Introduction	15
1.1. General	15
1.2. Therapy possibilities	15
1.3. Finally	15
Chapter 2 - Safety	15
2.1. Preliminary notes	15
2.2. General	15
2.3. Exclusions	15
2.4. Current density	15
2.5. Electromagnetic interference	16
2.6. Product liability	16
2.6.1. Limitations of liability	16
Chapter 3 - Indications and Contra-Indications ..	17
3.1. Indications	17
3.2. Therapy	17
3.3. Contra-indications	17
Chapter 4 - Installation	17
4.1. Connection	17
4.2. Connection of the mains adapter	17
4.3. Switching on and self test	17
4.4. Disconnecting from the mains	17
4.5. Installation	17
4.6. Electromagnetic interference	18
Chapter 5 - Operation	18
5.1. Controls	18
5.1.1. Device	18
5.2. Operation	19
5.2.1. Setting by hand of the Endomed 182 ..	19
5.2.2. Pre-programmed protocols selection ..	20
5.2.3. Programming	20
5.3. Actions in relation to treatment	21
5.4. Connection and disconnection reactions ..	21
5.5. Electrolytic effects	21
Chapter 6 - Maintenance	21
6.1. Technical maintenance	21
6.2. Cleaning of the apparatus	21
6.3. Cleaning of accessories	21
6.3.1. Electrodes and sponge pads	21
6.3.2. Patient cable	22
6.4. Environment information	22
Chapter 7 - Fault conditions	22
7.1. Display's fail to light up	22
7.2. Code ErrX on the display	22
7.5. Insufficient or no output current	22
Chapter 8 - Specifications	23
8.1. Specifications	23
8.2. Technical specifications	23
8.3. Classification	24
Chapter 9 - Order detail	24

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1.1. General

The Endomed 182 is a unit for electrotherapy. The available currents are well suited for pain management in both acute and chronic conditions. The currents are easy selectable and adjustable. During the treatment all parameters are continuously controllable and visible

The equipment is only meant to be used by competent personnel in physiotherapy, rehabilitation or adjacent areas.

1.2. Therapy possibilities

Electrotherapy with both classic and isoplanary interferential current. Bipolar interferential current and asymmetric biphasic pulsed current (TENS) with or without the use of a BURST frequency.

Through the use of these currents it is possible to diminish muscle tone and pain. Furthermore the vascularisation and regeneration of the tissue can be improved.

1.3. Finally

You have made a wise choice in selecting the Endomed 182. We are confident that your unit will continue to give satisfaction over many years of use. Nevertheless, if you have any queries or suggestions, please contact your authorised distributor.

2.1. Preliminary notes

It is important that you read these operating instructions carefully before using the Endomed 182. Please make sure that these instructions are available to all personnel who operates the equipment.

Pay attention to the following before using the Endomed 182:

1. Keep yourself informed of the contra-indications (see chapter 3).
2. Never leave the patient unattended during treatment.
3. The apparatus may not be used in close proximity (i.e. less than 2 metres) to shortwave equipment.
4. The apparatus may not be used in so-called "wet rooms" (hydrotherapy rooms).

The manufacturer cannot be held responsible for the results of using this apparatus for any purposes other than described in these operating instructions.

2.2. General

When the unit is switched on all vital functions are controlled via an internal microprocessor. During the functioning of the unit the actual output of the current is measured continuously and compared to the requested output. Whenever a faulty situation is detected the unit will switch off the current immediately.

2.3. Exclusions

Electrotherapy is not intended for the following specific applications:

- Intercranial treatment, like sleep therapy;
- Transcordial treatments (applying electrodes at the thorax may increase the risk of heart fibrillation);
- Cervical occipital treatments.

2.4. Current density

The particular standard for Therapeutic Electrical Stimulation Equipment is IEC 60601-2-10 (BS 5724 part 2.10). The maximum permissible Current Density is 2 mA R.M.S. per cm². To verify that this value is not exceeded in electrotherapy, divide the current output in mA by the effective electrode area in cm²; a value of less than two should result.

Note that placing the electrodes should be done carefully and solidly, and that the use of small electrodes in combination with high intensities may cause skin irritations or burns.

2.5. Electromagnetic interference

Simultaneous connection of a patient to High Frequency Therapy (shortwave or microwave) equipment and an electrical stimulator must not be done as this may result in burns at the site of the stimulator electrodes. Therefore, simultaneous connection is not allowed. For more information see chapter 4 Installation.

2.6. Product liability

A law on Product Liability has become effective in many countries. This Product Liability law implies, amongst other things, that once a period of 10 years has elapsed after a product has been brought into circulation, the manufacturer can no longer be held responsible for possible shortages of the product.

2.6.1. Limitations of liability

To the maximum extent permitted by applicable law, in no event will Enraf-Nonius or its suppliers or resellers be liable for any indirect, special, incidental or consequential damages arising from the use of or inability to use the product, including, without limitation, damages for loss of goodwill, work and productivity, computer failure or malfunction, or any and all other commercial damages or losses, even if advised of the possibility thereof, and regardless of the legal or equitable theory (contract, tort or otherwise) upon which the claim is based. In any case, Enraf-Nonius's entire liability under any provision of this agreement shall not exceed in the aggregate the sum of the fees paid for this product and fees for support of the product received by Enraf-Nonius under a separate support agreement (if any), with the exception of death or personal injury caused by the negligence of Enraf-Nonius to the extent applicable law prohibits the limitation of damages in such cases.

Enraf-Nonius can not be held liable for any consequence resulting from incorrect information provided by it's personnel, or errors incorporated in this manual and / or other accompanying documentation (including commercial documentation)

The opposing party (product's user or it's representative) shall disclaim Enraf-Nonius from all claims arising from third parties, whatever nature or whatever relationship to the opposing party.

USA only:

Safety Precautions

This unit operates with high voltages. Servicing of the Endomed should be referred to qualified service technicians or returned directly to the distributor. To assure continued compliance with FDA 21 C.F.R. 1050.10 standards, the Endomed should be calibrated and safety tested at least once each year. The producer recommends that all service work be performed by the distributor. However, a service manual can be obtained by writing to the factory for those who are qualified to service the device and are fully aware of their responsibilities.

Caution:

Use of controls or adjustments or performance of procedures other than those specified herein may result in hazardous exposure to ultrasonic energy.

See also the therapy book that can be ordered separately.

3.1. Indications

Diagnostics, electro palpation for:

- pain points;
- trigger points;
- hyperaesthetic regions;
- motor trigger points.

3.2. Therapy

A. Pain control in:

- pain points;
- trigger points;
- hyperaesthetic regions.

B Autonomic disorders such as:

- shoulder/hand syndrome;
- Raynaud's disease;
- Buerger's disease;
- Sudeck's dystrophy;
- neurological diseases;
- myalgia.

C. Muscle stimulation:

- for restoration of kinaesthesia;
- for treatment of atrophy;
- for stimulation of the internal and external sphincters in the treatment of urinary incontinence;
- revalidation;
- muscle strengthening (sport);
- neurological diseases
- spasm inhibit.

D. Iontophoresis.

3.3. Contra-indications

- pyrexia;
- tumours;
- tuberculosis;
- localised inflammation;
- thrombosis;
- pregnancy;
- pacemakers;
- metal implants.

4.1. Connection

- Mains supply connections must comply with the national requirements regarding medical rooms.
- Prior to connection of this apparatus to the mains supply, check that the voltage and frequency stated on the type plate correspond with the available mains supply.
- The adapter is a part of the supply circuit on which the device's safety partly depends. The approvals for the Endomed 182 are only valid if used in combination with this type of adapter ENB-1530.

It is not permitted to connect the Endomed 182 to another type of adapter than type **ENB-1530**.

4.2. Connection of the mains adapter

- Connect the supplied mains adapter to the connector [14].
- Connect the mains adapter to a wall socket. The flashing dot in the upper display [3]. indicates that the unit is connected to the mains supply and stand-by.

4.3. Switching on and self test

- Switch On the unit, using the On/Off button [1]
- Immediately after switching on, the unit carries out a self test.

Check that:

- All LED's light up briefly.
- At the end of the test a beep can be heard.

Contact your supplier if this is not the case.

4.4. Disconnecting from the mains

- Switch the Endomed 182 Off, using the button [1].
- Pull out the mains plug of the mains adapter from the wall socket.
- Disconnect the mains adapter from the Endomed 182.

4.5. Installation

- Do not install the unit in a location near to a heat source such as a radiator.
- Avoid exposure to direct sunlight, rain, excessive dust, dampness, mechanical vibrations and shocks.
- Cordless telephones should not be used in close proximity of the unit.
- This unit should not be used in so-called "wet-rooms" (hydrotherapy rooms).
- The unit has to be installed in such way that liquid cannot enter.

Should any liquid ingress into the housing of this equipment, unplug the unit from the wall socket (if connected) and have it checked by an authorised expert.

4.6. Electromagnetic interference

- Operation in close proximity (i.e. less than 2 metres) to working shortwave or microwave therapy units may produce instability in the Endomed 182 output.
- To prevent electromagnetic interference, we strongly advise using separate mains groups (phases) for the Endomed 182 and the shortwave or microwave equipment, and a distance of at least 2 metres is kept between the Endomed 182 and the shortwave or microwave equipment.
- Ensure that the mains cable of the shortwave/microwave unit does not come near the Endomed 182 or the patient.

If problems with electromagnetic interference persist, please contact your authorised distributor.

CHAPTER 5 - OPERATION

5.1. Controls

(See jacket flap)

5.1.1. Device

[1] On / off button

This button is to switch the unit on and off.

[2] Current form display

This display shows the current selection by LED. To change the current form push button [16].

[3] Display for AMF and modulation frequency (spectrum)

This display shows the readout from AMF or modulation frequency. The display mode is changeable with button [4]

[4] Selection button

With this button you change the read-out mode in display [3] from AMF to Modulation frequency

[5] Adjust buttons for AMF and Modulation frequency

Use this [-] / [+] buttons to adjust the AMF and Modulation frequency.

[6] Display for timer and indicators for modulation program

Read-outs display for timer and the indications LED's for the modulation programs.

[7] Selection button for modulation program

Selection button for the modulation programs. Selection will be indicate by LED in display [6]

[8] Adjust buttons for the treatment time

Use this [-] / [+] buttons to adjust the timer

[9] Display for intensity read-out in channel 1

Read-outs for intensity in channel 1, channel selection and balance.

[10] Display for intensity read-out in channel 2

Read-outs for intensity in channel 2

[11] Selection button channel 1 or 2

This button is to select the channel for adjusting the intensity with the buttons [13].

[12] Connector for patient cable.

Connection for the patient cable for both channels.

Connection of accessories other than ones specified by the manufacturer can adversely affect the safety of the patient and correct functioning of the equipment, and is therefore not permitted.

For combined applications use only Enraf-Nonius type BF equipment. The very low leakage current of this type of equipment ensures absolutely safe therapy.

[13] Adjusting buttons for the intensity and pre-programmed protocols

Use this [+] / [-] buttons to adjust intensity and program selector for the programs P0 until P30.

[14] Connector

Connect here the main adapter.

[15] Memory button

To activate the memory function, recall and store button. Read-out is in display [9] en [10].

[16] Selection button Current forms

To select the current forms. The selection will be displayed on display [2]

[17] Connection plug remote control

Connect here the optional remote control.

Connection of a remote control other than one specified by the manufacturer can adversely affect the safety of the patient and correct functioning of the equipment, and is therefore not permitted.

5.2. Operation

5.2.1. Setting by hand of the Endomed 182

Switch on

- Connect the standard main adapter to connector [14]
- Connect the main adapter plug to the mains socket.
- Switch the unit on with the button [1]

During the start-up period runs a self-test.

- Connect the patient cable adapter on connector [12]
- Fix the electrodes to the patient.

Therapy selection

- Select the current form with button [16]
- 2P: Pre-modulated interference
- 4P: Interference classical
- : Isoplanair interference
- : BI-asymmetrical biphasic pulsed current (TENS)
- : BURST tens (Bi-asymmetrical biphasic pulsed current)

The selected current will be indicated on display [2] by green LED

Adjusting AMF

- Select the AMF read-out with button [4]
- The LED light-up next to 
- Adjust the AMF frequency with the buttons [5]

Read-out is on Display [3]

Adjusting Modulation frequency

- Select modulation frequency read-out with button [3]
- The LED light-up next to 
- Adjust the modulation frequency with the buttons [5]

Read-out is on display [3]

Adjusting the modulation program

- Adjust the modulation program with button [7]. (3 settings are possible 1-1, 6-6, 1-30)

The selection is indicated in display [6] by LED next the figure 1-1, 6-6 or 1-30

Adjusting treatment time

Adjust the treatment time with the buttons [8]

Read-out is on display [6]

Adjusting the intensity in channel 1 and 2

- Select the channel with button [11]
- Selected channel will be indicate by LED in display [10]
- Adjust the intensity with the buttons [13]

Countdown of treatment timer starts.

Repeat this handling for channel 2.

*If interference classical or isoplanair is selected. The intensity will be adjusted automatically in both channels. Use the channel selection for balance.

Stop the treatment

- Adjust the timer to zero or select a other current form

5.2.2. Pre-programmed protocols selection

Enraf-Nonius stored parameters on the P1 until P30. The treatment suggestions are described on the therapy card.

Switch on

- Connect the standard main adapter to connector [14]
- Connect the main adapter plug to the mains socket.
- Switch the unit on with the button [1]

During the start-up period runs a self-test.

- Connect the patient cable adapter on connector [12]
- Fix the electrodes to the patient.

Program selection

- Push the button [15]
- Select a program number(P 1 - P 30) with the buttons [13] read-out in display [9 + 10]
- After the selection with the buttons [13] parameters will be displayed of the selected program.

Start the treatment

- The program will start after the adjusting of the intensity.
- The countdown starts.

Stop the treatment

Adjust the timer to zero or select a other current form

5.2.3. Programming

The program P0 can be programmed. This position is ideal to store the most frequent used parameters. When the unit start-up these parameters of program P0 will be displayed. The program number P1 – P30 can be overwritten at the same way.

Switch on

Connect the standard main adapter to connector [14]
Connect the main adapter plug to the mains socket.
Switch the unit on with the button [1]

Program selection

- Push the button [15].
- Select a program number(P 1 - P 30) with the buttons [13] read-out in display [9 + 10].
- After the selection with the buttons [13] parameters will be displayed of the selected program.

Treatment parameters: Adjusting or changing

- Adjust the parameters as decrypted under 9.1 Settings of the Endomed182” (follow the steps until adjusting the intensity).

Storing

- Push at the program button [15] during minimal 3 seconds.
- The parameter will be stored under the selected program number. The storing will be confirmed by a double beep.

5.3. Actions in relation to treatment

Before treatment

- Check the patient for absolute and relative contraindications.
- Test the heat sensibility of the treatment area.
- If the skin is hairy shaving is the best.

Connection of accessories other than ones specified by the manufacturer can adversely affect the safety of the patient and correct functioning of the equipment, and is therefore not permitted.

5.4. Connection and disconnection reactions

Constant Current (CC) output characteristics may cause unpleasant connection and disconnection reactions if the electrodes are not securely placed or lose contact with the skin. Make sure the output is set to 0 mA when you apply or remove the electrodes.

5.5. Electrolytic effects

Electrolysis occurs under the electrodes when current types with a DC component are applied. Because the greatest concentration of electrolytic by-products caused by ion migration occur under the electrodes, we recommend the use of the supplied sponges to keep the effects to a minimum. Make sure that the sponges are kept well moistened and place the thick side of the sponge between the flexible rubber electrodes and the patient.

6.1. Technical maintenance

We recommend to have the unit checked annually. This may be done by your supplier, or by another agency, authorised by the manufacturer. It is also recommended that a record of the service history is kept for all activities relating to service and maintenance. In some countries this is even obligatory.

Maintenance and all repairs should only be carried out by an authorised agency. The manufacturer will not be held responsible for the results of maintenance or repairs by unauthorised persons.

The service and/or maintenance must be carried out conform the procedure described in the service manual of the unit.

Opening of the equipment by unauthorised agencies is not allowed and will terminate any claim to warranty.

6.2. Cleaning of the apparatus

Switch off the unit and disconnect it from the mains supply. The apparatus can be cleaned with a damp cloth. Use lukewarm water and a non-abrasive liquid household cleaner (no abrasive, no alcohol content solution).

6.3. Cleaning of accessories

6.3.1. Electrodes and sponge pads

The rubber electrodes should be cleaned with lukewarm water. For stubborn stains or dirt, cleaning and disinfection of the rubber electrodes can be done using a 70% alcohol solution. During this process, the black colouration of the rubber electrodes may be stained. This does not affect the operation of the electrodes.

The sponge pads must be washed after use in warm water with a household cleaner. Then they must be rinsed with water, drained thoroughly and then dried, after they have been placed like roofing tiles. Damaged sponge pads must be replaced.

In areas where the water is soft, it is possible that the water may conduct poorly resulting in abnormally low currents. In this case it is possible to use a saline solution to improve the conductivity of the water.

When the electrodes are not being used, the sponge pads should be removed. This increases the operational life of the rubber electrodes. We advise that you should keep an extra set of sponge pads in stock.

6.3.2. Patient cable

Check the cable regularly for damages and/or bad electrical contact. We also advise that you keep an extra patient cable in stock.

6.4. Environment information

Your Endomed 182 contains materials which can be recycled and/or are noxious for the environment. At the end of the duration of life, specialised concerns can take apart the Endomed 182 and sort out the noxious materials and materials for recycling. By doing so you contribute to a better environment.

Please ensure that you are well informed of the local rules and regulations regarding to the removal of equipment and accessories.

CHAPTER 7 - FAULT CONDITIONS

7.1. Display's fail to light up

- Check if the mains adapter is connected to the unit and mains supply.

7.2. Code ErrX on the display

The unit has discovered a fault using its "self-test". Remove any applicators or cables from the output sockets from the equipment, switch the unit back on. If the code re-appears, contact your supplier. The unit is probably defective.

- * ErrX: number is mentioned on display [6], running from 1 until and including 9.

7.5. Insufficient or no output current

- Check the patient cable for breaks or poor contact.
- Ensure that the sponges are sufficiently damp.
- If necessary, use a salt solution.
- Clean the electrodes.
- Check for poor contact.

8.1. Specifications

Current

Intensity	: 0 - 99 mA
Timer	: 0 - 60 minutes

2P : 2 pole mid- frequency current :

Carrier wave	: 4 kHz
AMF frequency	: 5 - 200 Hz
Modulation frequency	: 0 - 180 Hz
Modulation program	: 1-1, 6-6, 1-30 seconds

4P : Interference classical and isoplanairy

Carrier wave	: 4 kHz
AMF frequency	: 5 – 200 Hz
Modulation frequency	: 0 – 180 Hz
Modulation program	: 1-1, 6-6, 1-30 seconds

Bi – asymmetrical biphasic pulsed current (TENS)

Pulse width	: 150 μ s
Frequency	: 1 – 200 Hz
Modulation frequency	: 0 – 180 Hz
Modulation program	: 1-1, 6-6, 1-30 seconds
Burst frequency	: 2 Hz

8.2. Technical specifications

Mains adapter

Type mains adapter	: ENB-1530
Mains voltage	: 100 - 240 Volt
Frequency	: 50/60 Hz.
DC Output	: 15 Volt
Max. Output current	: 2 A

Unit

Medical device classification	: IIb (according to Medical Device Directive (93/42/EEC)
Safety class	: II* type BF**, according to IEC 60601-1
Patient leakage current	: Better than IEC requirement ($IEC \leq 100 \mu A$)
Ditto, first fault condition	: Better than IEC requirement ($IEC: \leq 500 \mu A$)
Safety tests	: CE-MDD (TÜV Rheinland)
Dimensions	: 22 x 20 x 20 cm (w x d x h)
Weight	: 1.8 kg.

Environment conditions for transport and Storage

Environment temperature	: -10° till +70° C
Relative humidity	: 10 till 100 %
Atmospheric pressure	: 500 till 1060 hPa

Environment conditions normal use

Environment temperature	: 10° till 40° C
Relative humidity	: 10 till 90 %
Atmospheric pressure	: 500 till 1060 hPa

8.3. Classification

Medical class IIb

This equipment complies with all requirements of the Medical Device Directive (93/42/EEC).

**International safety norm**

The unit is confirming the international norm for medical electronic devices IEC 60601-1.

*II : Safety class II (double insulated).

**BF : The equipment has a floating patient circuit. The Endomed 182 meet the requirements of IEC 60601-1.

Technical modifications reserved

CHAPTER 9 - ORDER DETAIL

For the ordering data of the Endomed 182, standard accessories and additional accessories refer to the Catalogue for Physiotherapy and leaflet.

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